


TheraCell



TheraFuze DBF[®]

FIBER BOAT[™]



FIBER BOAT[™]

1 LEVEL • 2 LEVEL • 3 LEVEL BOATS

TheraFuze DBF[®]

Fiber form for Posterolateral Fusion with Fiberlok™ technology.

ONE LEVEL BOAT - ORDER: TC-FB-SM



TWO LEVEL BOAT - ORDER: TC-FB-LG



THREE LEVEL BOAT - ORDER: TC-FB-XL



What is it?	A demineralized bone fiber (DBF) 'boat' that can be loaded with autograft and placed in the posterolateral gutters for one, two or three level fusions.
What is it made of?	Demineralized bone fibers molded into a boat shape with FiberLok™ technology. 100% cortical bone, with no excipients. Produced in compliance with HCT/P regulations under Title 21 Code of Federal Regulations (CFR), Part 1271.
How can it be used?	The lyophilized boat may be rehydrated with saline or the patient's blood or bone marrow. The surgeon may then use the depression in the boat to hold autograft.
How does it work?	TheraCell's proprietary FiberLok™ technology provides strength to the hydrated boat to allow it to be manipulated and placed in the intended surgical site. Osteoconductive with osteoinductive potential to stimulate long term bone formation.
What sizes are available?	SM: One level ● LG: Two level ● XL: Three level
Patents	The implants and/or techniques associated with this product may be covered by one or more of the United States Patents: US 9,572,912 US 9,486,557 AUS 2014253753 Other patents pending

TO ORDER OR FOR MORE INFORMATION PLEASE CONTACT US AT - 630-953-9594

THERACELL, INC.

14930 VENTURA BLVD., SUITE #325
SHERMAN OAKS, CA 91403

630-953-9594

www.theracellinc.com

TheraFuze DBF® is a registered trademark of TheraCell, Inc.
Fiber Boat™ and FiberLok™ are trademarks of TheraCell, Inc.

ALL RIGHTS RESERVED.

